



**Manual for Multi-purpose
Measurement Software,
Beans**

- Please note that the contents of the manual are subject to change without notice.



Field Logic Inc.

CONTENTS

Chapter1: Structure of Beans

| | |
|----------------------------|---|
| 1-1 Home Screen..... | 4 |
| 1-2 Graph Screen | 4 |
| 1-3 Report Screen | 5 |
| 1-4 Monitor Screen | 5 |
| 1-5 Settings1 Screen | 6 |
| 1-6 Settings2 Screen | 6 |
| 1-7 Settings3 Screen | 7 |

Chapter2: Basic Usage

| | |
|------------------------------------|----|
| 2-1 START/STOP | 8 |
| 2-1-1 START..... | 8 |
| 2-1-2 STOP..... | 8 |
| 2-2 Screen Changeover | 9 |
| 2-3 Check Measurement Status | 9 |
| 2-4 CSV Save | 12 |
| 2-5 Report Print-outs | 12 |

Chapter3: Make Simple Changes in Basic Settings

| | |
|--|----|
| 3-1 Home | 13 |
| 3-1-1 Change a Channel Item on Home | 13 |
| 3-1-2 Change a Display Title | 14 |
| 3-1-3 Change Display Image..... | 14 |
| 3-1-4 Change a Data Format..... | 15 |
| 3-2-1 Relative Graph | 15 |
| 3-2 Graphs..... | 15 |
| 3-2-2 Absolute Graph | 17 |
| 3-3-1 Changeover to Daily/ Monthly/ Yearly Report..... | 21 |
| 3-3 Report..... | 21 |
| 3-3-3 Change/ Add/ Delete | 21 |
| 3-3-2 Change Date of Report | 21 |
| 3-3-4 Default Display of Report Item..... | 25 |

Chapter4: Change Measurement Settings


| | |
|--|-----------|
| 4-1 Lock/ Unlock Protection | 26 |
| 4-1-1 Unlock Protection | 26 |
| 4-1-2 Unlock Protection | 26 |
| 4-2 Settings1 Settings | 27 |
| 4-2-1 Change Sampling Interval/ Data Save Interval | 27 |
| 4-2-2 Measurement Channel Settings | 28 |
| 4-3-1 Add/ Delete Formula | 31 |
| 4-3 Settings2 | 31 |
| 4-3-2 Change Formula | 33 |
| 4-4 Settings3 | 34 |
| 4-4-1 Change Communication System Information | 34 |
| 4-4-2 Start-up Settings | 34 |
| 4-4-4 Start-up Settings | 35 |
| 4-4-3 Change Measuring Target | 35 |

Chapter5: Appendix

| | |
|---------------------------------|-----------|
| 5-1 Tables | 36 |
| 5-1-1 Settings1 Items | 36 |
| 5-1-3 Functions | 37 |
| 5-1-2 Arithmetic Operator | 37 |

Chapter1: Structure of Beans


1-1 Home Screen



Shows measurement items and values.

Such converted values as amount of CO2 reduced can be displayed as well.

1-2 Graph Screen



Shows graphs of items displayed on Home.

Possible to choose either a relative or an absolute graph.

1-3 Report Screen

The screenshot shows a software window titled "9Forest Beans". The menu bar includes Home, Graph, Report (which is selected), Monitor, Settings1, Settings2, and Settings3. The main area displays a table with the following data:

| | Irr. [kWh/m ²] | Temp. [°C] | DC Pwr [kWh] | AC_Pwr [kWh] | Amount of CO2... |
|-----|----------------------------|------------|--------------|--------------|------------------|
| 0' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 1' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 2' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 3' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 4' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 5' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 6' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 7' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 8' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 9' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 10' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 11' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 12' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 13' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 14' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 15' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 16' | 1.07 | 22.0 | 649.4 | 442.6 | 139.2 |
| 17' | 0.57 | 22.0 | 346.3 | 236.1 | 74.2 |
| 18' | 0.00 | 0.0 | 0.0 | 0.0 | 0.0 |
| 19' | 0.00 | 0.0 | 0.0 | 0.0 | 0.0 |
| 20' | 0.00 | 0.0 | 0.0 | 0.0 | 0.0 |
| 21' | 0.00 | 0.0 | 0.0 | 0.0 | 0.0 |
| 22' | 0.00 | 0.0 | 0.0 | 0.0 | 0.0 |
| 23' | 0.00 | 0.0 | 0.0 | 0.0 | 0.0 |
| Sum | 18.81 | 16.5 | 11305.8 | 7760.3 | 2440.6 |

Buttons at the bottom include AUTO-SET, DATE (D/M/Y), UPDATE, PRINT, CSV, RENEW, START, and STOP.

Designed for print-outs of reports, and save as CSV.
Report items can be added or deleted either manually or automatically.

1-4 Monitor Screen


The screenshot shows a software window titled "9Forest Beans". The menu bar includes Home, Graph, Report, Monitor (selected), Settings1, Settings2, and Settings3. The main area displays two tables:

| Real-time Channel Values | | | | Real-time Formula Values | | | | |
|--------------------------|-----------------|--------------|----------|--------------------------|-------------|--------------|-------|--------|
| No. | Name | Raw Data | Meas... | No. | Name | Formula | Value | Unit |
| 6 | AC_Voltage... | 420.000 | 420.000 | P0 | Reduced_CO2 | ch010*0.3145 | 2.296 | kg-CO2 |
| 7 | AC_Current(R) | 2.000 | 2.000 | P1 | Reduced_OIL | ch010*0.227 | 1.657 | l |
| 8 | AC_Current(S) | 2.000 | 2.000 | | | | | |
| 9 | AC_Current(T) | 2.000 | 2.000 | | | | | |
| 10 | AC_Power | 7.300 | 7.300 | | | | | |
| 11 | AC_Energy(...) | 1188.000 | 1188.000 | | | | | |
| 12 | AC_Energy(...) | 12345678.000 | 12345... | | | | | |
| 13 | Irradiance | 4.000 | 4.000 | | | | | |
| 14 | Temperature | 18.499 | 18.499 | | | | | |
| 15 | PV_Temp | 4.196 | 4.196 | | | | | |
| 16 | Wind_Speed | 20.000 | 20.000 | | | | | |
| 17 | spare | 4.000 | 4.000 | | | | | |
| 18 | INV_Run | 1.000 | 1.000 | | | | | |
| 19 | INV_Wait | 1.000 | 1.000 | | | | | |
| 20 | INV_Stop | 1.000 | 1.000 | | | | | |
| 21 | INV_Trouble | 1.000 | 1.000 | | | | | |
| 22 | INV_GridTra... | 1.000 | 1.000 | | | | | |
| 23 | MesCount | 1.000 | 1.000 | | | | | |
| 24 | Latitude(N/S) | 0.000 | 0.000 | | | | | |
| 25 | Latitude(XX") | 35.000 | 35.000 | | | | | |
| 26 | Latitude(YY") | 46.824 | 46.824 | | | | | |
| 27 | Longitude(E/... | 0.000 | 0.000 | | | | | |
| 28 | Longitude(XX... | 129.000 | 129.000 | | | | | |
| 29 | Longitude(Y... | 48.782 | 48.782 | | | | | |
| 30 | Vessel_Speed | 20.200 | 20.200 | | | | | |
| 31 | BowDirection | 120.000 | 120.000 | | | | | |

Buttons at the bottom include RUN, RENEW, START, and STOP.


Shows all real-time values of both measured and calculated.

1-5 Settings1 Screen




Shows settings of all measurement items.

1-6 Settings2 Screen



Shows settings of formulas.

1-7 Settings3 Screen



Shows network-related settings, measuring target setting, and password settings.


Chapter2: Basic Usage

2-1 START/STOP

2-1-1 START

Click on START button. (①)

Once measurement starts, the indicator on the left bottom turns on as it is shown on the right.
(②)



2-1-2 STOP



Click on STOP. (①)

A confirmation dialog will pop up, and click OK. (②)

Type in password. (③)



*Password is required only when launched from startup.



2-2 Screen Changeover

Click on one of the tabs to change over to the other screens.

※ It will switch back to Home when there is more than a hour of no mouse control.



2-3 Check Measurement Status

Indicators:

Shows measurement status.

Different status are shown as follows:



: Green



Status : Normal



: Yellow



Status : Partial Failure



: Red



Status : Failure



: Black / No light on



Status : Measurement Off

*Examples of Possible Failure

*Error signals from measurement instrument

*Connection error with measurement instrument

Home:

Such Measured values as Irradiance and Alternative Current Power are shown on Home.


Check if all measurements are taken, and none of the values is either too high or low.

On the right is when measurement is normal.

※ Such measurement values as DC/AC Generated Power vary depending on rated output of PV system.

Examples of real-time values

*Displayed contents may vary depending on each setting.



Monitor:

Shows all real-time measurements as well as operation values.

Designed for morning measurements real-time.

If measurements are not normally displayed, there is possibility for a glitch in system.

Real-time values

| No | Name | Raw Data | Meas. |
|----|-----------------|---------------|---------------|
| 6 | AC_Voltage | 420.000 | 420.000 |
| 7 | AC_Current(R) | 2.000 | 2.000 |
| 8 | AC_Current(D) | 2.000 | 2.000 |
| 9 | AC_Current(T) | 2.000 | 2.000 |
| 10 | AC_Power | 7.300 | 7.300 |
| 11 | AC_Energy | 123456789.000 | 123456789.000 |
| 12 | AC_Energyf | 123456789.000 | 123456789.000 |
| 13 | Irradiance | 4.000 | 4.000 |
| 14 | Temperature | 10.499 | 10.499 |
| 15 | PV_Temp | 4.196 | 4.196 |
| 16 | Wind_Speed | 20.000 | 20.000 |
| 17 | spare | 4.000 | 4.000 |
| 18 | INV_Run | 1.000 | 1.000 |
| 19 | INV_Wat | 1.000 | 1.000 |
| 20 | INV_Stop | 1.000 | 1.000 |
| 21 | INV_Trouble | 1.000 | 1.000 |
| 22 | System_Start | 1.000 | 1.000 |
| 23 | MeCount | 1.000 | 1.000 |
| 24 | Latitude(NS) | 0.000 | 0.000 |
| 25 | Latitude(NS)_ | 36.000 | 36.000 |
| 26 | Latitude(YY)_ | 45.924 | 45.924 |
| 27 | Longitude(E)_ | 0.000 | 0.000 |
| 28 | Longitude(XX)_ | 129.000 | 129.000 |
| 29 | Longitude(YY)_ | 46.02 | 46.02 |
| 30 | Vertical_Signed | 20.200 | 20.200 |
| 31 | BowDirection | 120.000 | 120.000 |

| No | Name | Formula | Value | Unit |
|----|-------------|--------------|-------|--------|
| P0 | Reduced_CO2 | c#01010.3145 | 2.296 | kg-CO2 |
| P1 | Reduced_OIL | c#01010.227 | 1.657 | l |

Possibility in measurement fault

| No | Name | Raw Data | Meas. |
|----|---------------|----------|---------|
| 6 | AC_Voltage | 420.000 | 420.000 |
| 7 | AC_Current(R) | 0.000 | 0.000 |
| 8 | AC_Current(D) | 0.000 | 0.000 |
| 9 | AC_Current(T) | 0.000 | 0.000 |
| 10 | AC_Power | 0.000 | 0.000 |
| 11 | AC_Energy | 0.000 | 0.000 |
| 12 | AC_Energyf | 0.000 | 0.000 |
| 13 | Irradiance | 4.000 | 4.000 |
| 14 | Temperature | 10.499 | 10.499 |
| 15 | PV_Temp | 4.196 | 4.196 |
| 16 | Wind_Speed | 20.000 | 20.000 |
| 17 | spare | 4.000 | 4.000 |
| 18 | INV_Run | 1.000 | 1.000 |



Graphs:

Visually shows real-time measurements in graphs.

Users can choose to show measurements either in relative or absolute graph.

※ If graphs are not normally displayed, there is possibility for a glitch.

(e.g. There is an abnormally rapid fall(s) in data graphs; Max./ Min. values are both zero.)



2-4 CSV Save

Click on CSV on Report screen.

(①)


CSV Output dialog will pop up.

Select Measurement Items on report to save data for the selected report items or all measurement items. Select time format from Y for yearly, M for monthly and D for daily(②).


Choose save period from Term Selection(③).

Click Save (④), and choose where to save the file.

※In default setting, a file name contains a date on which the file is created. After naming it,click Save.



CSV output window




2-5 Report Print-outs

Click on PRINT. ()

Select a printer and execute it.

※Set a printer before using the print-out function.




Chapter3: Make Simple Changes in Basic Settings

3-1 Home

3-1-1 Change a Channel Item on Home

Double-click on one of the item boxes on Home. (①)

Home Tab Display Settings dialog will pop up.





Click on CH and select one of items from drop-down listbox.

(②)

Click OK.(③)

*When the change is made to CH, so is to the Title. However, the image will not change.


Home Tab Display Setting



3-1-2 Change a Display Title

Double-click on one of the item boxes on Home.
Home Tab Display Settings dialog will pop up.

Click on input area of Title, and type in a name(①).
Click OK to end.(②)



3-1-3 Change Display Image

Double-click on one of the item boxes on Home.
Home Tab Display Settings dialog will pop up.




Click on Image(①), and select one of images from drop-down listbox.

Click OK.(②)

As a result, the original image has been replaced with the selected image.

(Refer to the right)

*Clicking on “status” box(③) enables you to show an image when a chosen item is above or below zero.
Without checking the box, the image will not change in accordance with the value.




3-1-4 Change a Data Format


Click on Data Format radio button on Home Tab Display Settings (①), and choose one of the formats.

Click OK.(②)

*As for channels such as Fault ch, data format is fixed, and not to be changed.



Radio buttons




3-2 Graphs

3-2-1 Relative Graph


3-2-1-1 On/Off Relative Graphs

By ticking on / off the check boxes (①), it is possible to turn graphs on / off.




3-2-1-2 Changeover to Daily,Monthly, Yearly Graph

Click on one of the radio buttons to change over to the other type of graphs(①).



3-2-1-3 Change Date on Graph for Checking Past and Present Data

Click on an arrow (①) to change D for date/ M for month/Y for year. Date can also be changed by clicking <>>.




3-2-1-4 Relative Graph Display Applied/ Unapplied

Tick/ untick the Relative checkbox to apply/ unapply the display method. (①)

※ Ticked: Shows relative display graphs, using maximum and minimum values of data on a certain date. Suitable for tracking moves in data values.

※ Unticked: Shows relative display graphs, using data collected a month previously and after from start date. Suitable for tracking every-day moves in data values.




3-2-2 Absolute Graph

3-2-2-1 Absolute Graph Settings

Click on SETTINGS (①)

to show Graph Settings window.





On absolute graph, it is possible to show graphs up to four.

On the setting window, click on tabs to change from one axis to another.

Change axis title and unit, if necessary. Set axis title and unit for a real-time value where it says "Real-time", and another title and unit for accumulative value where it says "Others".

Make settings on graph channels in Legend area.

Click on a CH list box (①), and select a channel to be displayed from the drop-down list.




If more than one channel has already been set, only channels, of which unit is same as the existing channel, are displayed on the list.

※One the right is where DCW has been set, and the only item left to choose is ACW.

Click an item either on Legends (Real-time) or Legends (Others) (①), and change graph description(s).


Set a legend for a real-time value where it says "Real-time", and another one for a cumulative value where it says "Others".



Set max. and min. value for each graph.

※If not set at all, both values will automatically be set.


Click on OK (①) to reflect the settings.



Clicking on AUTO-SET clears all settings that have been made, and reflect default settings.

3-2-2-2 Changeover to Daily , Monthly , Yearly Graph

Click on the list box (①), and choose Real-time, D (Daily), M (Monthly), or Y (Yearly).




3-2-2-3 Change Date on Graph for Checking Past and Present Data

Click on an arrow (①) to change D for date/ M for month/ Y for year.


Date can also be changed by clicking<<>>.

To display the latest data, click on RENEW (②).



3-2-2-4 Print-Out

Click on PRINT (①).



Select your printer (②) and print graphs (③).

※ Make sure that at last a printer is connected to your PC.


※ Relative graphs can not be printed.



3-3 Report

3-3-1 Changeover to Daily/ Monthly/ Yearly Report

Click one of three radio buttons to choose from daily/ monthly/ yearly report. (①)




3-3-2 Change Date of Report

Select date from drop-down listbox (①), or click right/ left arrow button. (②)



3-3-3 Change/ Add/ Delete

To Change an Item: Double-click on an item. (①)
(In this case, temperature is chosen.)



Report Item Settings dialog will pop up.

Select an item from Channel No.

(②)

In this case, Irr. is selected.

Click OK. (③)



Check that an item selected is on Report.

Click on RENEW button if the change was made correctly.




An Item changed

| Int. [kWh/m ²] | Tmp. [°C] | DC Pwr [kWh] | AC Pwr [kWh] | Amount of CO2... |
|----------------------------|-----------|--------------|--------------|------------------|
| 0° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 1° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 2° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 3° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 4° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 5° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 6° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 7° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 8° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 9° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 10° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 11° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 12° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 13° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 14° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 15° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 16° 1.0° | 22.0 | 49.4 | 442.6 | 139.2 |
| 16° 0.5° | 22.0 | 46.3 | 236.1 | 74.2 |
| 18° 0.00 | 0.0 | 0 | 0 | 0 |
| 19° 0.00 | 0.0 | 0 | 0 | 0 |
| 20° 0.00 | 0.0 | 0 | 0 | 0 |
| 21° 0.00 | 0.0 | 0 | 0 | 0 |
| 22° 0.00 | 0.0 | 0 | 0 | 0 |
| 23° 0.00 | 0.0 | 0 | 0 | 0 |
| Sum | 18.81 | 16.5 | 1995.8 | 2440.6 |

To Add an Item:

Right-click on * mark (④) and select 'Insert'. (⑤)



Report Item Settings dialog will pop up.


Click on Channel No. and select an item on the list. (⑥)
e.g. In this case, IRR is chosen.

Finally, click on RENEW to reflect the settings.

Report Item Setting




An Item added



Right-click on an existing item (⑦) to insert a new item into records. The new item will be inserted next to the existing item just selected.


Select "Insert". (⑧)

Choose one of channel numbers , and click on RENEW button.



To Delete an Item:

Right-click on an item (⑨), and select "Erase". (⑩)



Finally click RENEW button.

3-3-4 Default Display of Report Item

Click on AUTO-SET button.(①)




Confirmation message will pop up, and click on OK button.

(②)

Finally, click on RENEW button.

Confirmation Window



Chapter4: Change Measurement Settings

4-1 Lock/ Unlock Protection

4-1-1 Unlock Protection



Click on STOP button to terminate measurement.

※ Not possible to unlock protection while measuring.

Double-click on the key mark

(①), and type in password (②), and click OK. (③)

Once it is done, the key will be open. (Refer to the right)



4-1-2 Unlock Protection

Double-click on the open key mark.(①)



The key will be closed.

(Refer to the right)



4-2 Settings1 Settings

4-2-1 Change Sampling Interval/ Data Save Interval

※ In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)



Click on an edit box of sampling interval/ data save interval (①), change the existing number in the box (①), and click on RENEW button. (②)
The number of measurements will change automatically.

※ Data save interval has to be multiples of sampling interval.
If not, an error message will pop up, and settings will be reloaded.

※ Click on NO button to load previously saved settings (④).

※ The number of measurement is automatically set.
※ Data save interval needs to be the integral multiple of sampling interval.
※ Default Settings:
※ Sampling Interval; 6 secs
※ Data Save Interval; 60 secs
※ The smaller the data save interval, the lower the performance of data display and calculation it may become.

Sampling interval,data save interval and the number of measurements



4-2-2 Measurement Channel Settings

Unlock protection.



Mode:

Double-click on a mode (①),
and select a mode from
drop-down listbox.

Name:


Double-click on a name, and type
a new name in. (②)

※The new name will be reflected
on the same names on the
other tabs.



Saving Format:


Double-click on a saving format item (③), and select a format from drop-down listbox.



Unit:

Double-click on a unit item (④), and select a unit from drop-down listbox.

If there is no unit to choose in the listbox, type a new unit in. The new unit will automatically be saved.




Mask:

Tick/ untick a mask checkbox.

(⑤)


※If max. and min. values have to be set, you need to tick it.



Max. and Min. Values:

Double-click on Max/ Min Value of masked item, and type a new value in. (⑥)

※Mainly used for analogue measurements.




Identification Value:

Double-click on an ID Value item, and select a new item from dropdown listbox. (⑦)

Renew:


Click on RENEW button (⑧) and OK button. (⑨)



Undo Settings:

Click on NO button to undo settings.

(⑩)



4-3 Settings2



4-3-1 Add/ Delete Formula

Add a Formula Pattern1:

※ In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)

To insert a formula, or add a new one, double-click on the name item of *mark (①), and type a new name in (②).



※Refer to Appendix for available arithmetic operators and functions.



Next, double-click on Measurements item, (③) and type a new formula in. Double-click unit and identification value, and select one from drop-down listbox. (④)

Click on RENEW button. (⑤)


Click on OK. (⑥)




Add a Formula Pattern 2:

To insert a new item in between existing items, put the pointer on an item, onto which to insert a new one.

Then, right-click (⑧) and select Insert. (⑨)




The selected item will be copied and inserted.



Delete an Item:

To delete an item, put the pointer onto an item to delete, (⑩) right-click and select Erase. (⑪)

Click RENEW.



4-3-2 Change Formula

※ In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)

Unlock protection.

(Refer to 4-1-1)


Double-click where to make changes at. (①)

If it is Name, double-click on Name and type a new name in.

If it is Measurements, double-click on Measurements and type a new formula in.

If it is Unit, double-click on Unit and select one from drop-down listbox.

If it is Identification Value, doubleclick on IDV and select one from dropdown listbox.




Finally, click on RENEW button.

(②)

Message of "Settings saved" will pop up, and click OK. (③)

※Click NO button to undo previously saved settings. (④)

※It may take a few seconds to save settings.



4-4 Settings3



4-4-1 Change Communication System Information

※ In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)

There are settings for Site ID/ Database Communication Port/ Data-base Server.

※ Changes to the settings are not recommended.


After making changes, click RENEW button (①). Message of "Settings saved" will pop up, and click OK. (②)



4-4-2 Start-up Settings

※ In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)

Tick/ untick "Minimize and display on the task tray" checkbox. ()
Click RENEW button.
Lock protection.
(Refer to 4-1-2)





4-4-3 ChangeMeasuring Target

※In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)

Select a target from General-purpose,Photovoltaic Generation, or Photothermal Utilization. (①)

Click RENEW button. (②)

Message of "Settings saved" will pop up, and click OK (③) or click CANCEL (④) to go back to previously saved settings.



4-4-4 Start-up Settings

※ In order to change settings, protection needs to be unlocked.
(Rfr: 4-1-1 Unlock Protection)

Click on Password Setting box and type a new one in. (①)

Click RENEW button (②).



Chapter5: Appendix

5-1 Tables

5-1-1 Settings1 Items

| <Item> | <Subitem> | <Overview> |
|----------------------|-------------------|--|
| Saving Format | Integrated | Obtains integrated value out of data within Data Save Interval time. |
| | Avg. | Obtains average value out of data within Data Save Interval time. |
| | Nearest | Obtains nearest value out of data within Data Save Interval time. |
| Unit | A | Direct/ Current Ampere |
| | V | Direct/ Current Voltage |
| | KW | Direct/ Current Power |
| | KW/m ² | Irradience |
| | Kg-CO2 | CO2-Coverted Value |
| | °C | Temperature |
| | Times | The Number of Measuremunts |
| Mask | Mask | Switch on/ off setting for max./ min. value. |
| Identification Value | TIME_a | The Number of Measurements |
| | IRR_h | Irradiance |
| | TEMP_m | Temperature(Avr.) |
| | DC_V_mS | Direct Current Voltage(Avr.) |
| | DC_A_m | Direct Current Ampere(Avr.) |
| | DC_W_h | Direct Current Power |
| | AC_V_m | Alternative Current Voltage(Avr.) |
| | AC_A_m | Alternative Current Ampere(Avr.) |
| | AC_W_h | Alternative Current Power |
| | CO2_h | CO2 Reduced |
| | _h | Quantity(AVr. Value per hour) |
| | _m | Avr. |
| | _a | Integrated |
| | _b | Status |

5-1-2 Arithmetic Operator

| <Operator> | <Overview> | <E.G.> | <Answer> |
|------------|----------------------------|---------|----------|
| + | Sum | 2+3 | 5 |
| - | Difference | 2-3 | -1 |
| * | Multiplication | 2*3 | 6 |
| / | Division | 4/2 | 2 |
| % | Residue | 5%4 | 1 |
| ^ | Power | 2.0^3.0 | 8 |
| l/ | Square Root | l/25.0 | 5 |
| ll/ | Cube Root | ll/27.0 | 3 |
| ! | Factorial | 5! | 120 |
| !! | Factorial(Prefix Operator) | !!5 | 120 |
| @ | Absolute Figure | @-5.0 | 5 |
| & | Binary AND | 91 & 15 | 11 |
| | Binary OR | 31 3 | 35 |
| # | Binary XOR(Exclusive OR) | 17 # 5 | 20 |
| ~ | Binary NOT | ~1 | -2 |
| << | Binary Left Shift | 1<<4 | 16 |
| >> | Binary Right Shift | 8>>2 | 2 |

5-1-3 Functions

| <Function> | <Overview> | <E.G.> | <Answer> |
|------------|---|------------------|--------------------|
| abs(z) | Absolute Value | abs(-17.4) | 17.4 |
| cbrt(d) | Cube Root | cbrt(27) | 3 |
| ceil(x) | Roundup(The smallest Integer larger than argument) | ceil(42.2) | 43 |
| degrees(d) | Radian against d Degree | degrees(0.5) | 28.6478897565412 |
| exp(x) | Exponent | exp(1.0) | 2.7182818284590452 |
| floor(x) | Cut-down(The largest integer smaller than argument) | floor(42.2) | 42 |
| in(x) | Natural Logarithm | in(2.0) | 0.6931471805599453 |
| log(x) | Natural Logarithm larger than 10 | log(100.0) | 2.000000000000000 |
| log(y,x) | Logarithm of X larger than y | log(2.0,64.0) | 6.000000000000000 |
| mod(l,j) | Residue of i/j | mod(9.2,4) | 1.2 |
| pi() | Circle Ratio | pi() | 3.14159265358979 |
| pow(y,x) | y to the x Power | pow(2,4) | 16 |
| radians(d) | Degree against x Radian | radians(45) | 0.785398163397448 |
| random() | Random Digits between 0.0_1.0 | | |
| round(x) | Round-off | round(42.2) | 42 |
| round(x,i) | Round off to i | round(42.2765,2) | 42.28 |
| setseed(d) | Set Seed of Random Digit | | |
| sign(x) | Code for Argument x (-1,0,1)) | sign(-8.4) | -1 |
| sqr(x) | Square Root | sqr(2.0) | 1.414213562373095 |
| trunc(x) | Truncation | trunc(42.2) | 42 |
| trunx(x,i) | Round down to i decimal place | trunx(42.2765,2) | 42.27 |